

# RYERSON

DEPARTMENT OF CHEMISTRY, BIOLOGY AND CHEMICAL ENGINEERING  
FACULTY OF ENGINEERING AND APPLIED SCIENCE

Mr. Ravi Kanipayor  
Aquasure Technologies  
1259 Ernest Ave.  
London, ON


November 30, 2000

Dear Mr. Kanipayor,

This letter is to confirm that I have processed coliform/*E.coli* samples with both your STEPPI incubator and a traditional incubator. Using the Colilert-18 kit (IDEXX) and an inoculum of 1000CFU/ml, I was able to read the test results after 17-18 hours as outlined by the manufacture. The results were as reliable as those obtained with a standard incubator. Furthermore, I found that your incubator maintained a constant temperature of 35°C through out the experiment.

If you require further information, please contact me.

Sincerely,



Kim Gilbride, PhD